




Creating balance in healthcare



Amfora Shower stretcher



Amfora Shower stretcher



6300 4010 / 6200 4015

Amfora Shower stretcher hydraulic height adjustable

- 70 x 160 cm.
- Synthetic collection tray with drain hose.
- 4 stainless steel double lock swivel castors.
- Adjustable backrest.
- Side rails.
- Maximum weight 135 kg.

Also available in:

- 70 x 190 cm.

6200 4011 / 6200 4016

Amfora Shower stretcher electric height adjustable

- 70 x 160 cm.
- Synthetic collection tray with drain hose.
- 4 stainless steel double lock swivel castors.
- Adjustable backrest.
- Side rails.
- Maximum weight 135 kg.

Also available in:

- 70 x 190 cm.





6300 4030 / 6300 4035
**Amfora Shower stretcher
wall mounted, hinged**

- 70 x 160 cm.
- Folding legs.
- Pole fixation.
- Adjustable backrest.
- Side rail.
- Maximum weight 135 kg.

Also available in:

- 70 x 190 cm.



6300 4040 / 6300 4045
**Amfora Shower stretcher
wall mounted, hinged**

- 70 x 160 cm.
- Synthetic collection tray with drain hose.
- Folding legs.
- Pole fixation.
- Adjustable backrest.
- Side rail.
- Maximum weight 135 kg.

Also available in:

- 70 x 190 cm.

6200 4050 / 6200 4055
**Amfora Shower stretcher
fixed height**

- 70 x 160 cm.
- Synthetic collection tray with drain hose.
- 4 stainless steel double lock swivel castors.
- Adjustable backrest.
- Side rails.
- Maximum weight 135 kg.

Also available in:

- 70 x 190 cm.



6200 5005 / 6200 5010

Amfora Shower Stretcher
electric height adjustable

- 70 x 160 cm.
- Wall mounted.
- Synthetic collection tray with drain hose.
- Foldable, gas spring assisted.
- Lowest position, 43 cm.
- Adjustment range of 43-103 cm.
- Suitable for bathrooms with a minimum height of 230 cm.
- Synthetic screen guard.
- Adjustable backrest.
- Side rail.
- Maximum weight 135 kg.
- **Optional:** extra battery pack.

Also available in:

- 70 x 190 cm.



For more information, visit our website: www.lopital.com

Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk, Postbus 56, 5060 AB Oisterwijk
Tel +31 (0)13 5239300, Fax +31 (0)13 5239301, E-mail info@lopital.nl, www.lopital.nl

