



Creating balance in healthcare



# Luxal Shower stretcher



6100 2200

# Luxal Shower stretcher

6100 2200

## Luxal Shower stretcher

This predominantly aluminium frame shower stretcher, electrically height adjustable, offers extra comfort to both the resident / client and caregiver with its comfortable form, simple to operate hand switch and shock free height adjustment.

Access to both sides of the stretcher enables the caregiver to move freely around the stretcher. The significant reduction in weight makes the stretcher very easy to move and suitable to maneuver around tight spaces.



*Height adjustment 55 - 95 cm.*



*Tilting position for water drainage.*



- Electric adjustable height.
- Simple to operate hand switch.
- Includes battery charger (wall mounted).
- Battery is easy to change.
- Optional extra battery.
- Emergency stop.
- Emergency lowering device.
- Includes mattress and pillow.
- Tiltable to drain the water.
- Folding down safety rails.
- Central lockable swiveling castors with one directional castor.
- Clearance space up to 20 cm.
- Swivel castors 15 cm.
- **Optional:** extra battery pack.





1. Directional wheel.



2. Neutral / drivable position.



3. Central complete lock position.



Manual control & battery.



Discharge hose.

## Power / Operation

High / low drive:	24VDC-lifting pin, max. 5A
Manual control:	4-key handset with hook, IP66
Control box:	24VDC, max. 250W
Emergency brake:	push button with locks and emergency stop switch
Battery:	24V, 2,9Ah, with quick-changing chuck
Battery charger:	230V AC, electrical power 100W, wall mounted
Emergency descent switch:	stick shift, directly from the control box
Complete electrical system:	IP65
Wheels:	diameter 15 cm
Side rails:	turnable, requires both hands when opening

## Dimension / weight

Overall size:	204 x 78 x 75/115 cm (L x W x H)
Top surface area:	190 x 65 cm
Mattress height:	20 cm
Minimum height:	55 cm
Maximum height:	95 cm
Water discharge hose:	inwendig Ø4 cm
Wheel base:	141 cm
Wheel gauge:	65 cm
Clearance space:	20 cm
Weight (nett):	67 kg
Weight capacity:	180 kg maximum

For more information, visit our website: [www.lopital.com](http://www.lopital.com)

Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk, Postbus 56, 5060 AB Oisterwijk  
 Tel +31 (0)13 5239300, Fax +31 (0)13 5239301, E-mail [info@lopital.nl](mailto:info@lopital.nl), [www.lopital.nl](http://www.lopital.nl)

