



**ChildReflex**



**51005610  
ChildReflex**

Caring for a disabled child while using a traditional shower-toilet chair, is very cumbersome. The caregiver will be forced to stand in an unnatural position and has to cope with static back strain. ChildReflex is the total solution to help the disabled child shower and with toilet needs. ChildReflex provides comfort and safety to the caregiver as well as the child, thanks to the stepless electrical tilting device which provides optimal safety. When set at the highest level, caregiver and child are able to make eye contact.



*Optional: Neoprene inlay to accommodate small/medium/maximum sitting area.*

**There are three different ergonomically designed showerseats available that offers a comfortable sitting area to the child.**

- Small sitting area. Specifications: 24 cm wide, 32 cm deep. **51005610**  
Age group: 3-6 years old.
- Medium sitting area. Specifications: 28 cm wide, 34 cm deep. **51005611**  
Age group: 7-10/11 years old.
- Large sitting area. Specifications: 34 cm wide, 36 cm deep. **51005612**  
Age group: 10-14/15 years old.

(These are approximate age guidelines)