

# AUVS Disinfection Box

KR615L

## AUVS Disinfection Box

The UV Box - uses Photon Trap Technology to disinfect smaller handheld and various devices in less than a minute! The UV Box protects staff, patients and their families.





Place any item inside The UV Box without stacking or overlapping



Close the door and turn knob to lock



Press the red button to begin the disinfection cycle



### Characteristics

- Proven efficacy
- FDA Consulting Lab Tested
- Environmentally friendly (Reduce cleaners and wipes)
- Durable
- Economical
- Mobile and accessible
- Bulb lasts 16,000 hours or 50+ years under normal use

Overall Dimensions	10" H x 10" D x 16" W
Enclosure Dimensions	3.5" H x 8" D x 15.5" W
Weight	10 lbs.
Power	110V, AC

Hospitals invest \$100,000 or more in UV robots, foggers and an unending supply of wipes, but spend little to no money addressing the worst culprits of cross contamination: The mobile devices that move freely in, out and around your hospital.

*Cornerstone Medical Inc. brings you the industry's best products from companies who focus their R&D to meet healthcare's evolving needs.*



Cornerstone Medical Inc., 4-889 Barton Street, Stoney Creek, ON L8E 5V1  
 T 905-945-2522 or 1-800-652-3895 | F 905-945-5611 | E info@cornerstone-medical.com

[www.cornerstone-medical.com](http://www.cornerstone-medical.com)