



**CORNERSTONE
MEDICAL**
Innovative solutions | Client-focused service

Discover with **StoryPanels™**

THE SILENTIA SCREEN SYSTEM





"Planets with illustrations for older children", article no. 6109-APO

*StoryPanels in treatment rooms and wards
places medical equipment in the background
and brings children to more playful and
positive thoughts.*





Unfold the magic of StoryPanels

A picture is worth a thousand words and could even have a healing impact. Silentia's StoryPanels offer a new way to give children a playful and positive atmosphere in the healthcare sector by creating a positive distraction.

A child in a hospital is far from home and among unfamiliar faces. But then a nurse starts talking about the colorful picture with planets and funny figures. The child starts using its fantasy and imagination. Changing focus to something new and positive.

Art has a stress-reducing effect by allowing patients to think about something other than their condition. StoryPanels takes Silentia's folding screens to a new level by its combination of creating both privacy and a more positive atmosphere. A powerful tool for the healthcare sector to comfort children.





Design that stimulates the imagination


"Each detail in the design is based on research and experience. The colors, shapes and funny figures all have a purpose in order to give a positive psychological effect, a calm feeling and to stimulate the imagination," says Stefanie Schöffmann, Austrian designer who has developed the StoryPanels' design for Silentia. The result is a playful and inspiring design which helps caregivers to make children feel more comfortable.

Stefanie Schöffmann • Designer, CEO look! design



The "Planets only" design, here in combination with panels in the color Ocean, article no. 6107-APP-O





StoryPanels offer health-oriented spaces for children and create a less clinical look.



A new solution to improve patients' well-being

The healthcare sector is more and more aware of the benefits of using art and images in health care. Art allows patients to focus on something other than their condition. Studies have shown that pictures can influence well-being and relieve patient pain and stress.

Nurses and caregivers can use the *StoryPanels* as an easy tool ready at hand to distract, stimulate and play during treatment sessions. It comforts and encourages the children and parents and thereby provides a better working environment. The design blends in easily with the interior without distracting nurses and staff, who see it every day.

Children can look at the motifs from many positions – lying down, sitting up, close by or from a distance. They will feel more calm and comfort and by looking at the motifs again and again they will discover new details, new stories and new fantasies. Also older children can explore the textures and clever details.

Parents sitting beside their child find comfort in exploring the textures, getting lost in details and daydream some of their worries away. Or make up new, little bedside stories to encourage their child.


The last 10 years Silentia has worked together with designers and nurses on how to use folding screens printed with pictures to improve patients' well-being. With *StoryPanels* Silentia has developed a print technology and standardized system to offer health-oriented spaces for children and create a less clinical look.

*StoryPanels is a new way to create a
playful environment for the child,
parent and healthcare professional
to relieve a stressful situation.*




Choose between three planet designs or make a mix


Planets only




Planets with illustrations for younger children, 1–5 years old




Planets with illustrations for older children, 6 years and up




Panel height




Choose your Silentia Screen solution. See the whole range at silentia.us




Mobile Bed End Screen



Mobile Folding Screen



Wall mounted



Floor mounted

Find height and number of panels that fits your Silentia Screen solution *

FOLDING SCREENS NO. OF PANELS	ART. NO. PLANETS ONLY	ART. NO. YOUNGER CHILDREN	ART. NO. OLDER CHILDREN	LENGTH	HEIGHT	WEIGHT
5	6105-APP(-X)	6105-APY(-X)	6105-APO(-X)	4'1"	4'9"	18.7lbs
	6205-APP(-X)	6205-APY(-X)	6205-APO(-X)	4'1"	5'5"	20.9lbs
	6305-APP(-X)	6305-APY(-X)	6305-APO(-X)	4'1"	6'1"	25.4lbs
7	6107-APP(-X)	6107-APY(-X)	6107-APO(-X)	5'9"	4'9"	26.5lbs
	6207-APP(-X)	6207-APY(-X)	6207-APO(-X)	5'9"	5'5"	30.9lbs
	6307-APP(-X)	6307-APY(-X)	6307-APO(-X)	5'9"	6'1"	35.3lbs
9	6109-APP(-X)	6109-APY(-X)	6109-APO(-X)	7'5"	4'9"	34.2lbs
	6209-APP(-X)	6209-APY(-X)	6209-APO(-X)	7'5"	5'5"	39.7lbs
	6309-APP(-X)	6309-APY(-X)	6309-APO(-X)	7'5"	6'1"	45.2lbs
11	6111-APP(-X)	6111-APY(-X)	6111-APO(-X)	9'	4'9"	40.8lbs
	6211-APP(-X)	6211-APY(-X)	6211-APO(-X)	9'	5'5"	47.4lbs
	6311-APP(-X)	6311-APY(-X)	6311-APO(-X)	9'	6'1"	54lbs
13 ¹⁾	6113-APP-X	6113-APY-X	6113-APO-X	10'8"	4'9"	48.5lbs
	6213-APP-X	6213-APY-X	6213-APO-X	10'8"	5'5"	56.2lbs
	6313-APP-X	6313-APY-X	6313-APO-X	10'8"	6'1"	64lbs
15 ¹⁾	6115-APP-X	6115-APY-X	6115-APO-X	12'4"	4'9"	56.2lbs
	6215-APP-X	6215-APY-X	6215-APO-X	12'4"	5'5"	65.0lbs
	6315-APP-X	6315-APY-X	6315-APO-X	12'4"	6'1"	74lbs
BED END SCREENS	6522-APP	6522-APY	6522-APO	4'11"	5'1"	52.9lbs
	6523-APP	6523-APY	6523-APO	4'11"	6'1"	59.5 lbs
	6525-APP	6525-APY	6525-APO	6'7"	5'1"	62.8 lbs
	6526-APP	6526-APY	6526-APO	6'7"	6'1"	71.6 lbs

The table shows our standard sizes in 20" length increments. Screens are also available in increments of 10" if the total length is between 29.5" and 12'3". The standard heights – 4'9", 5'5" and 6'1" – can be lowered by – 3.9" by determining the height of the legs of the wheels. The height is decided when you order your Silentia Folding Screen.

Note: The longest screens (12–15 sections) are intended for use in an L-shape around the bed.

Important: With a X in the end of the Art No. the screen is combined with laminate panels according to Daylight Duo (DD) configuration and X represents the color symbol for laminate White, Ocean, Forest, Savannah, Tropical, Birch, Lime, Mandarin or Daylight.

¹⁾ Available only as Daylight Duo, see page 7.

* Ask Silentia or your local distributor for advice

Easy to combine

Combine privacy with a positive atmosphere. StoryPanels come alive due to the light effect that shines through the Silentia Daylight panel.

StoryPanels offer full flexibility and scalability with endless variations. They come in six original panels suitable for three different heights. The design works with either Daylight screen panels or color panels.

Choose between three layer options or make a mix:

- Planets only (APP)
- Planets with illustrations for younger children (APY)
- Planets with illustrations for older children (APO)

StoryPanels™ is

- made with a hard, smooth and easy to clean surface
- made to withstand all disinfectants
- durable and UV-resistant

EasyClean™ for optimal hygiene

The screens have smooth surfaces that withstand all disinfectants used in care situations today. The enclosed wheels can easily be removed and washed separately in a dishwasher. All products meet strict hygiene requirements and all surfaces are accessible for cleaning – that's EasyClean.

EasyReturn™ folds back easily with one hand

The patented hinges have a built-in memory, which enables easy handling and prevention of cross infection

using only one contact point: the ergonomically designed handle.

EasyClick™ – fixed or mobile, the choice is yours

Movability is one click away. A wall-mounted screen can be mounted on a trolley for mobile use, or to any other of our system's components. There are many possibilities, making this a complete system, easy to use in all departments.


See all the many possibilities at silentia.us



"Planets with illustrations for older children" with Bed End Screen, article no. 6522-APO

Combine StoryPanels with color panels.*

Daylight Duo screens must be used when the screen exceeds eleven panels.



*Article no. 6111-APP-O
Example of a combination with the color Ocean*

Article no. 6109-APP-O

Article no. 6107-APP-O

** Ask Silentia or your local distributor for advice*



Discover more
about **StoryPanels™**
and the entire range of
The Silentia Screen System
at silentia.ca



Cornerstone Medical Inc. | 10-889 Barton Street | Stoney Creek, ON L8E 5V1
Phone: 905.945.2522 or toll-free 1.800.652.3895 Fax: 905.945.5611
info@cornerstone-medical.com | www.cornerstone-medical.com



THE FUTURE IN PRIVACY & HYGIENE SOLUTIONS

Silentia AB • Box 108, SE-311 22 Falkenberg, Sweden • Tel. +46 (0)346-485 80
Fax +46 (0)346-485 89 • info@silentia.se • www.silentia.se



EN ISO 9001:2015



EN ISO 14001:2015

MADE IN SWEDEN

Our products are copyright-protected in various ways, including registered designs, patents and patents pending.
We reserve the right to make changes without prior notice. For more information, please visit our website.