

Tango XL Shower Commode Chair

5100 5705

Tango XL Shower and Commode Chair

The Tango XL Shower and Commode chair provides patient care and bathing in a responsible way. The Tango XL provides optimum comfort for the patient as well as the caregiver.






Infection Control Friendly
 Designed for Easy Cleaning
 Smooth Edges and Surfaces



Characteristics

- Smooth and Easy Drivability
- Designed for easy cleaning
- Lockable plastic swivel castors, stainless steel wheel housings
- PUR seat with toilet opening
- Rail for toilet bucket or bedpan
- Removable and swingaway footrests
- Can be wheeled over the toilet
- Rilsan coated frame
- 250kg/551lb
- Optional: Frontlocking armrests

Dimensions (cm)



Cornerstone Medical Inc. brings you the industry's best products from companies who focus their R&D to meet healthcare's evolving needs.



Cornerstone Medical Inc., 10-889 Barton Street, Stoney Creek, ON L8E 5V1
 T 905-945-2522 or 1-800-652-3895 | F 905-945-5611 | E info@cornerstone-medical.com

www.cornerstone-medical.com